

Product Release Notes

Product: PONEMAH Physiology Platform
Model: P3 Plus
Version: 4.60-SP2 (Service Pack 2 for version 4.60)
Build: J03263 (CD Build)
Date: October, 2007

Product Release Notes for PONEMAH Physiology Platform (P3 Plus) Version 4.60-SP2 indicate revisions made to the P3 Plus core application since release of P3 Plus Version 4.60-SP1. For information regarding changes to the software from previous versions, please refer to the Release Notes folder located on the Version 4.60 CD. Product Release Notes indicate only revisions to application contents that are part of CD Part #J03263 – Build Version 4.60-SP2.

Product Release Notes for P3 Plus do not include information regarding revisions of the P3 Plus Analysis Modules.

Notice for organizations that must comply with FDA's Good Laboratory Practices (GLP) and 21 CFR Part 11 Electronic Records; Electronic Signatures: P3 Plus Versions may contain **Preview Features**. These **Preview Features** are listed in the Product Release Notes table under the column, "Type of Change". A **Preview Feature** indicates that enhancements have been made to P3 Plus, but have not been validated. Instead, Data Sciences International (DSI) has opted to delay complete validation until receiving comments from customers regarding use of these features. Further validation of these features will be performed in later releases of P3 Plus. There may be additional **Preview Features** that had been documented in previously released versions that are not documented here. These features are not available unless manually enabled by the user. If documentation is needed regarding these features, please contact the Technical Support Group at DSI.

Reference #	Type of Change	<p>Key: N = New Feature; E = Enhancement; F = Fix</p> <p>Description</p>
Study Protocol Option		
2740	F	<p>The Study Protocol Option is a feature that provides structured handling and naming of data collected. Backup and Restore features allow set up information and run information to be archived, backed up, or restored on workstations.</p> <p>Restoring setup information could cause an issue if a study is created and data is collected, and then archived and deleted, and the original setup information restored to continue future collections. This could result in incorrect file names and error messages.</p> <p>Error checking has been added to ensure proper functionality if setup information must be archived, deleted and routinely restored during a study.</p>
2598	F	<p>Reviewing data through study could cause Ponemah to close when saving a marks information. If attributes are set such that many derived/calculated parameters do not report data ("x" listed in column), saving a marks section could result in the Ponemah program closing. This has been corrected to account for absent, calculated values.</p>
Data Manager		
2741	F	<p>When deleting files from the hard disk using Data Manager, file paths containing 85 characters or greater were not properly handled and could have resulted in the Ponemah program closing. This has been corrected to handle the Window maximum path length of 260 characters.</p>