

Workflow for programming infusion protocol, downloading, activation for RCV Variable Infusion Protocol.

RCV Variable

BookMarks of UserManual

iPRECIO Global Workflow Pg. 9 & 10
Review Checklist for precautions on Pg. 8


Starting a new study:
with number of groups/animals/drug concentration etc.
Pg. 30

Infusion & Flow Rate modes Review
Pg. 25 - 27
Study Programming
Pg. 31 - 35

Pump Detection (Recognition) & Program Pump
Pg. 36 & 37 respectively

Surgical Guide, Initial Fill & Activation
Pg. 40, 43 & 44 respectively

iPRECIO Management
Pg. 61




Header Step

Start a New Study Go to Header Step of iPRECIO® Software. Fill in details. All fields marked with *(asterisk is required) and these are:

- Required Fields**
1. Study Name
 2. Compound Name or ID
 3. Number of Groups
 4. Number of Animal in the group
 5. Study Number/ID
 6. Study Date
 7. Operator
 8. Concentration

- Example**
1. RCV Variable Workflow
 2. Ang II
 3. 2
 4. 5
 5. RCV Variable Workflow Example
 6. 10/08/03
 7. John Smith (Select your name)
 8. 100

Comment
For this example 2 groups, one for saline and second for Ang II. 5 animals per group. It is also recommended to fill all fields as it will serve as a study record.




Groups Step

Animal Groups Setting Need to fill in the following fields:

- Required Fields**
1. Concentration (if not filled in Header Step)
 2. Estimated Max/Min Animal Weight
 3. Enter Group Number/ID
 4. Select Infusion Mode
 5. Flow Rate Mode
 6. Select Dose or Flow Rate


- Example** (Group2, Ang II)
1. 100
 2. 275 max & 225 min
 3. Saline Group & Ang II Group
 4. RCV
 5. Variable
 6. Flow Rate

Comment
RCV variable requires a few more parameters than Instant Constant as an Initial Saline infusion is programmed for recovery period (RCV) prior to start of drug administration. Group 1 same as Instant Constant Setup. Not shown here.



Animals Step

Required fields to input are described here after. (Select the Ang II group.) Complete Animals Setting, then set Modes and Detection.



Animal Setting, Pump Assignment (Pump Detection/Recognition)

- Required Fields**
1. Animal ID, 2. Weight, 3. Sex
 4. Mode (Click **Set** button. [Post Recovery Mode])
 - a. Start of Administration Time/Date
 - b. Remaining tube length of iPRECIO pump
 - c. Additional catheter (if attached) ID and length or Measured Dead Volume
 - d. Flushing rate of Dead Volume,
 - e. RCV Start Time/Date,
 - f. RCV Flow rate
 - (NEXT to continue to Mode [Variable Flow Mode])
 - g. Step 1 flow rate and duration,
 - h. Step 2 flow rate and duration,
 - i. Step n flow rate and duration... (n=10 is max) Click on <Enter> once completed.
 5. Detect (Click **Detect** button to assign pump.)

* Repeat steps 1-5 for each animal/pump. Alternatively use the copy function to copy from the infusion protocol from the row above and then modify animal weight/ID as required. Once all pumps have been assigned to each animal, protocols may be downloaded to the pumps.

Example [Ang II Group with 3 steps used] 1. AG1, 2. 225, 3. F

Post Recovery Mode

Animal ID: AG1

Start of Administration (yy/mm/dd hh:mm): 12/12/12 09:00 ← 4.a.

Dead Volume Setup

Remain Tube Length: 60.0 (mm) ← 4.b.

Fundamental Dead Volume: 10.820 (ul)

Choose reference value (Dead Volume)

Measurement value / Calculated value

Inside Diameter: 0.0 (mm) Length: 0.0 (mm) ← 4.c.

Actual Dead Volume of additional catheter: 4.937 (ul)

Total Dead Volume: 30.004 (ul)

Flushing rate of Dead Volume: 30.0 (ul/hr) ← 4.d.

RCV Setup

RCV Start (yy/mm/dd hh:mm): 12/12/05 09:00 ← 4.e.

RCV Duration: 6 (days) 22.28

RCV Flow Rate: 3.0 (ul/hr) ← 4.f.

Refill Date / Time (if necessary): 12/12/12 07:29

Scheduled Exchange Date / Time (to): 12/12/12 07:58


Variable Flow Mode

Animal ID: AG5 Weight: 225(g)

Step	Flow Rate (ul/hr)	Duration (Day(s))	Repeat ON	Program No.	Times(n)	Start	End
Step1	5.0	7.0	← 4.g.			12/12/12 09:00	12/12/19 09:00
Step2	10.0	7.0	← 4.h.			12/12/19 09:00	12/12/26 09:00
Step3	15.0	7.0	← 4.i.			12/12/26 09:00	13/01/02 09:00
Step4	0.0	0.0					
Step5	0.0	0.0					
Step6	0.0	0.0					
Step7	0.0	0.0					
Step8	0.0	0.0					
Step9	0.0	0.0					
Step10	0.0	0.0					

Remaining Battery life: 5.30 (day(s))


*Repeat for all animals/pumps.




Pumps Step

Program Pump/Download, Pump Activation, Save Protocol and Monitoring


After all pumps are assigned to each animal, click **NEXT>>** to go Pumps window. Click **Program Pump** button.



Click the **Yes** button to download.




Download protocol to all pumps of all groups.



After finishing programming (downloading) to all group pumps, Status Manager of Pump is available.

To activate pumps, click **Activate** button in the Status Manager of Pump window. All pumps will need to be activated.

i Each pump must be activated at the appropriate <Activation Time Window>, 3 mins to 3 hours before infusion start time.



i Fill reservoir completely and ensure that solution reaches the distal end of iPRECIO pump catheter prior to activation.