

Species	Size	Weight (g)	Multicolor Kit PN	Single color kit PN (Qty 4)	Single collar PN	Color	Pitch (inches)	Thickness (inches)	Chamber PN
Mouse	Standard	20-25	601-6119-001	601-2153-001	009687-001	Blue	0.41	0.115	601-2026-001
				601-2154-001	009688-001	Green	0.38		
				601-2155-001	009689-001	Red	0.35		
				601-2156-001	011248-001	Yellow	0.32		
Rat	Standard	175-250	601-6106-001	601-2157-001	007685-001	Yellow	0.75	0.235	601-1307-001 or 010359-001
				601-2158-001	007684-001	Blue	0.70		
				601-2159-001	007683-001	Green	0.65		
				601-2160-001	007685-001	Red	0.60		
	Standard/Medium	250-325	601-6120-001	601-2161-001	007152-001	Yellow	0.85	0.235	601-3128-001
				601-2162-001	007151-001	Blue	0.80		
				601-2163-001	007150-001	Green	0.75		
				601-2164-001	007149-001	Red	0.70		
	Medium	325-400	601-6121-001	601-2165-001	007187-001	Yellow	0.90	0.235	601-0005-006
				601-2166-001	007186-001	Blue	0.85		
				601-2167-001	007185-001	Green	0.80		
				601-2168-001	007184-001	Red	0.75		
	Large	400-600	601-6122-001		009697-001	Yellow	1.00	0.235	601-2040-001
					009698-001	Blue	1.05		
					601-2117-001	Red	0.95		
	GP	Standard	260-320	601-6113-001		007816-003	White	1.02	0.235
320-380					007816-002	Orange	1.08		
380-440					007816-005	Red	1.15		
440-500					007816-001	Green	1.20		
500-560					007816-004	Blue	1.25		
560-620					011942-001	Grey	1.32		
620-680					011943-001	Yellow	1.39		
680-740					011944-001	Black	1.46		
Ferret	Standard	-	601-6123-001	601-6124-001	007816-002	Orange	1.08	0.235	601-2038-001
				601-6125-001	007816-004	Blue	1.25		
				601-6126-001	007816-001	Green	1.20		
				601-6127-001	007816-005	Red	1.15		

The pitch corresponds to this measurement :

