

30 September 2019

FinePointe Biosystem XA Study Import Capability

FinePointe V2.3.1.9 is the last version to offer validated support for the import of Biosystem XA data. Technical support and maintenance development for Biosystem XA was discontinued in January 2014. The removal of Biosystem XA data import support was done in consideration of this support discontinuation and the absence of communicated demand by customers.

Any attempts to conduct Biosystem XA data import with newer versions of FinePointe is done at the end user's own risk. Please note that compatibility issues may also exist with additional required system software, such as Microsoft® Windows®, MS SQL Server® and MS Office®.

It should also be noted that FinePointe V2.3.1.9 is considered a "Discontinued" version based on the DSI Software Support Lifecycle policy. Please the following link for additional details:

<https://www.datasci.com/products/software/dsi-software-support-lifecycle>

If Biosystem XA data must be accessible for an extended period of time due to government regulations and/or institutional requirements, it is recommended to maintain a full functional system, including computer, data acquisition hardware, and software, in archive to account for both software and computer hardware obsolescence.

If you have any questions, please contact DSI Technical Support at support@datasci.com or your local DSI support representative.



Michael W. Nielsen
Director, Scientific Services
Data Sciences Intentional